



# Skin Lesion Prediction:


A DERMATOLOGICAL COMPUTER-AIDED CLASSIFICATION PROJECT

IBM ADVANCED DATA SCIENCE SPECIALIZATION CAPSTONE

CHRISTOPHER LIEDEL, MAY 2021

# Technology


- ❑ IBM Watson Cloud and Jupyter notebooks.
- ❑ Python with data science libraries: pandas, scikit-learn, seaborn, matplotlib, numpy and keras.



# Data Used

**Dataset: HAM\_10000**  
(hmnist\_28\_28\_RGB.csv) for  
pigmented skin lesion  
diagnosis.

- ❑ Contains 10,000 images of skin lesions, both benign and malignant.
- ❑ A labels column.




# Data Used

## HAM\_10000\_metadata.csv file

Dataset has 7 possible diagnoses:

- ❑ 0: nv - Melanocytic nevi
- ❑ 1: mel - Melanoma
- ❑ 2: bkl - Benign keratosis-like lesions
- ❑ 3: bcc - Basal cell carcinoma
- ❑ 4: akiec - Actinic keratoses
- ❑ 5: vasc - Vascular lesions
- ❑ 6: df - Dermatofibroma


Also includes image/lesion ID, dx\_type, sex, age and location of lesion.



# Use Case: Multi-Class Classification

## Skin cancer is the most common human malignancy

- ❑ **Early detection** of skin cancer allows for a more successful and simpler cure.
- ❑ **Initial Diagnosis** is performed in a clinical setting limiting access to medical care for some.
- ❑ **Automating** skin lesion classification can be used as diagnostic aid in a clinical setting.
- ❑ **Automation** of skin lesion classification through mobile devices can potentially expand needed low-cost medical care outside the clinic.




# Data Analysis – diagnosis

## Lesions

**There are 7 types of lesions**


- ❑ 0: nv - Melanocytic nevi
- ❑ 1: mel - Melanoma
- ❑ 2: bkl - Benign keratosis-like lesions
- ❑ 3: bcc - Basal cell carcinoma
- ❑ 4: akiec - Actinic keratoses
- ❑ 5: vasc - Vascular lesions
- ❑ 6: df - Dermatofibroma

## Distribution






# Data Analysis – age & frequency



Histogram with fixed size bins (bins=18)

Value	Count	Frequency (%)
45	1299	13.0%
50	1187	11.9%
55	1009	10.1%
40	985	9.8%
60	803	8.0%
70	756	7.5%
35	753	7.5%
65	731	7.3%
75	618	6.2%
30	464	4.6%
Other values (8)	1353	13.5%

# Data Analysis – age/diagnosis/location




Value	Count	Frequency (%)
extremity	3195	24.2%
back	2192	16.6%
lower	2077	15.7%
trunk	1404	10.6%
upper	1118	8.5%
abdomen	1022	7.7%
face	745	5.6%
chest	407	3.1%
foot	319	2.4%
unknown	234	1.8%
Other values (6)	497	3.8%




# Data Analysis – sex/diagnosis/location

## Location Distribution by Sex



## Diagnosis and Sex




# Data Preparation

- ❑ Work on copies of the dataset (keep original intact)
- ❑ **Data Cleaning:**
  - Fill in missing values (e.g., mean, median, zero...) or drop the rows or columns.
  - Fix or remove outliers (optional).
- ❑ **Feature Selection** (optional):
  - Drop attributes that provide no usefulness for the task.
  - Drop duplicates.
- ❑ **Feature Engineering** where appropriate:
  - Discretize features
  - Decompose features (e.g., categorical).
  - Add transformations of features.
- ❑ **Feature Scaling:**
  - Standardize or normalize features


# Feature Analysis

## Between Variable Correlation Analysis

Phi k




Cramer's V




## Within Variable Distribution to target Correlation

Histograms and Scatter Plots




# Feature Analysis

**Using Mutual Information Scores to  
examine relationships between features**



**Class Label Distribution**




# Feature Transformation:

## Imbalanced Data


**Class Distribution**

```
Class=2, n=1099 (10.974%)  
Class=4, n=6705 (66.950%)  
Class=3, n=115 (1.148%)  
Class=6, n=1113 (11.113%)  
Class=5, n=142 (1.418%)  
Class=1, n=514 (5.132%)  
Class=0, n=327 (3.265%)
```



**Result of Using Synthetic Minority  
Oversampling Technique: SMOTE**

```
Class=2, n=6705 (14.286%)  
Class=4, n=6705 (14.286%)  
Class=3, n=6705 (14.286%)  
Class=6, n=6705 (14.286%)  
Class=5, n=6705 (14.286%)  
Class=1, n=6705 (14.286%)  
Class=0, n=6705 (14.286%)
```



# Performance Metric

We have a multi-class classification project so **accuracy** will be chosen for our performance metric, along with **precision**, **recall** and **specificity**.

**Macro** precision and recall used for balanced dataset.

**Weighted** precision and recall used for imbalanced dataset.

**90/10 Split** to Train and Test sets. Models evaluated on test sets.

**Model Evaluation** performed using both balanced and imbalanced data.

**Precision** — Out of all the examples that predicted as positive, how many are positive?

**Recall** — Out of all the positive examples, how many are predicted as positive?

## Used in Medicine

**Specificity** — Out of all the people that do not have the disease, how many got negative result?

**Sensitivity** — Out of all the people that have the disease, how many got positive test results?



# Tested Algorithms

Algorithm	Scores With Imbalanced Data	Scores With Balanced Data
Dummy Classifier	Acc .49    Weighted: Prec .49 / Rec .49 Specificity 0.861	Acc .15    Macro: Prec .15 / Rec .15 Specificity 0.858
Logistic Regression	Acc .69    Weighted: Prec .63 / Rec .69 Specificity 0.898	Acc .64    Macro: Prec .64 / Rec .64 Specificity 0.941
Decision Tree	Acc .63    Weighted: Prec .63 / Rec .62 Specificity 0.908	Acc .85    Macro: Prec .85 / Rec .85 Specificity 0.974
Random Forest Classifier	Acc .68    Weighted: Prec .62 / Rec .67 Specificity 0.892	Acc .92    Macro: Prec .92 / Rec .93 Specificity 0.986
XGBoost Classifier	Acc .73    Weighted: Prec .68 / Rec .73 Specificity 0.909	Acc .95    Macro: Prec .95 / Rec .95 Specificity 0.991
Neural Network with one hidden layer	Acc .73    Weighted: Prec .71 / Rec .73 Specificity 0.925	Acc .971    Macro: Prec .97 / Rec .97 Specificity 0.995
Neural Net with additional hidden layer of 64 neurons and dropout at ( 0.20)	Acc .73    Weighted: Prec .72 / Rec .74 Specificity 0.921	Acc .979    Macro: Prec .98 / Rec .98 Specificity 0.996
Neural Net with additional hidden layer of 64 neurons, dropout at ( 0.20) Tuned using Grid Search	NA	Acc .977    Macro: Prec .98 / Rec .98 Specificity 0.996

# Best Model: Neural Network

## Summary

Tuned, Fully Connected **Neural Network**

**Accuracy, Precision and Recall of 0.98**  
**Specificity of 0.996**

Epochs: 100

Batch size: 850

1st hidden layer: 512 neurons

2nd hidden layer: 64 neurons

1st activation: **relu**

2nd activation: **relu**

Dropout: 0.2


Output layer activation: **softmax**

Optimizer: **adam**

Loss: **sparse\_categorical\_crossentropy**


## Classification Report

	precision	recall	f1-score	support
0	0.99	1.00	0.99	665
1	0.99	1.00	0.99	676
2	0.95	0.98	0.97	670
3	1.00	1.00	1.00	621
4	0.97	0.91	0.94	695
5	1.00	1.00	1.00	705
6	0.95	0.98	0.96	662
accuracy			0.98	4694
macro avg	0.98	0.98	0.98	4694
weighted avg	0.98	0.98	0.98	4694




# Performance Charts

## Train and Test Performance



## Probability Curves



# Example of a prediction made using our best model

- ❑ The image used was selected from an external source and not taken from the dataset used to train and evaluate the model
- ❑ Model predicted actinic keratoses : class 4. This is the label given from the source as well.

```
from keras.preprocessing import image  
  
img_path = ('./Desktop/acct2.jpg')  
img = image.load_img(img_path, target_size=(158,158))  
img
```



```
# Classify image not in the dataset  
import cv2  
img = cv2.imread('./Desktop/acct2.jpg',0)  
  
img = cv2.resize(img, (58,58))  
cv2.waitKey(0)  
cv2.destroyAllWindows()  
  
print(class_label)  
  
prediction_cplx = tuned_model.predict_classes(img)  
  
print('The predicted class of the lesion is:', prediction_cplx[1])
```

```
{0: 'nv: melanocytic nevi', 1: 'mel: melanoma', 2: 'bkl: benign keratosis-like lesion', 3: 'bcc: basal_cell_carcinoma', 4:  
'akiec: actinic_keratosis', 5: 'vasc: vascular lesion', 6: 'df: dermatofibroma'}  
The predicted class of the lesion is: 4
```

# XGBoost

## Summary

Second best performing model is **XGBoost** with **accuracy**, **precision** and **recall** of **0.95**; **Specificity** of **0.991**


- ❑ max\_depth : 7
- ❑ booster: "gbtree"
- ❑ num\_classes: 7
- ❑ eval\_metric: "mlogloss"
- ❑ objective: "multi:softprob"

## Classification Report


xgb:	precision	recall	f1-score	support
0	0.96	0.99	0.98	665
1	0.96	0.98	0.97	676
2	0.92	0.92	0.92	670
3	0.99	1.00	1.00	621
4	0.95	0.79	0.86	695
5	1.00	1.00	1.00	705
6	0.86	0.94	0.90	662
accuracy			0.95	4694
macro avg	0.95	0.95	0.95	4694
weighted avg	0.95	0.95	0.95	4694

# XGBoost Performance

## Prediction Results



## Probability Curves





# Using the model

Trained model can be used in an application placed on a tablet device or smart-phone.

- ❑ A digital image is taken of the lesion in question.
- ❑ Image enters the transformation pipeline.
- ❑ The prepared image is fed into the pre-trained model.
- ❑ Model predicts the lesion class.
- ❑ The model prediction can be used as an aid to diagnosis.

Can be used outside of a clinical setting expanding access to needed low-cost medical care.

# Resources

Original Dataset Link:

<https://dataverse.harvard.edu/dataset.xhtml?persistentId=doi:10.7910/DVN/DBW86T>

Data Sources Link::

<https://www.kaggle.com/kmader/skin-cancer-mnist-ham10000>

**Image details** (Slide 5): Copyright: kali9 | Creator: kali9 | Royalty Free Source

**Jupyter Notebooks on IBM Cloud**

